

1722533012-08-V9



Part Number: 1722533012-08-V9

Manufacturer: Affinity Medical Technologies - a Molex company

Description: 8" PRE-CRIMP 1857/19 VIOLET

Data sheet:  1722533yyy Drawing

 Custom Jumper Part Numbering Guide

 Jumper Wires Length Guide

RoHs Status: Lead free / RoHS Compliant



Ship From: Hong Kong

Shipment Way: DHL/Fedex/TNT/UPS/EMS

Images are for reference only.  
See Product Specifications for product details.  
Buy 1722533012-08-V9 with confidence from  
Component-World.HK, 1 Year Warranty

Request For Quotation

| PRODUCT PARAMETER                                                                                                                                                                                                                                                                                                                                                                                                                                                             |                                                                         |                                  |                                                                                                                                       |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------|----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Part Number                                                                                                                                                                                                                                                                                                                                                                                                                                                                   | 1722533012-08-V9                                                        | Manufacturer                     | Affinity Medical Technologies - a Molex company                                                                                       |
| Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                   | 8" PRE-CRIMP 1857/19 VIOLET                                             | Lead Free Status / RoHS Status   | Lead free / RoHS Compliant                                                                                                            |
| Quantity Available                                                                                                                                                                                                                                                                                                                                                                                                                                                            | 75446 pcs                                                               | Data sheet                       | <a href="#">1722533yyy Drawing</a><br><a href="#">Custom Jumper Part Numbering Guide</a><br><a href="#">Jumper Wires Length Guide</a> |
| Category                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | <a href="#">Cable Assemblies</a>                                        | Wire Gauge                       | 18 AWG                                                                                                                                |
| Series                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | Ultra-Fit                                                               | Number of Rows                   | 1                                                                                                                                     |
| Number of Conductors                                                                                                                                                                                                                                                                                                                                                                                                                                                          | 1                                                                       | Moisture Sensitivity Level (MSL) | 1 (Unlimited)                                                                                                                         |
| Length                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | 8.00" (203.20mm)                                                        | Lead Free Status / RoHS Status   | Lead free / RoHS Compliant                                                                                                            |
| Detailed Description                                                                                                                                                                                                                                                                                                                                                                                                                                                          | Violet 18 AWG Jumper Lead Socket to Cable (Round) Gold 8.00" (203.20mm) | Contact Finish Thickness         | 30.0µin (0.76µm)                                                                                                                      |
| Contact Finish                                                                                                                                                                                                                                                                                                                                                                                                                                                                | Gold                                                                    | Contact End                      | Socket to Cable (Round)                                                                                                               |
| Color                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | Violet                                                                  |                                  |                                                                                                                                       |
| Component-World.com is a Reliable Stocking Distributor of Electronic Components. We specialize in all Affinity Medical Technologies - a Molex company series electronic components. We have 75446 pieces of Affinity Medical Technologies - a Molex company 1722533012-08-V9 in stock available now. Request a quote from electronics components distributor at Component-World.com, our sales team will contact you within 24 hours.<br>RFQ Email: info@Components-World.com |                                                                         |                                  |                                                                                                                                       |

| RELATED PRODUCTS                                                                   |                                                                                       |                                                                                |                |
|------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|--------------------------------------------------------------------------------|----------------|
|  | Part#: <a href="#">1722533012-08-Y7-D</a><br>Description: 8" PRE-CRIMP 1858/19 YELLOW | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-V7-D</a><br>Description: 8" PRE-CRIMP 1858/19 VIOLET | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-Y9-D</a><br>Description: 8" PRE-CRIMP 1857/19 YELLOW | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-S7</a><br>Description: 8" PRE-CRIMP 1858/19 SLATE    | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-W9</a><br>Description: 8" PRE-CRIMP 1857/19 WHITE    | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-S7-D</a><br>Description: 8" PRE-CRIMP 1858/19 SLATE  | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-R9</a><br>Description: 8" PRE-CRIMP 1857/19 RED      | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-V9-D</a><br>Description: 8" PRE-CRIMP 1857/19 VIOLET | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-R9-D</a><br>Description: 8" PRE-CRIMP 1857/19 RED    | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-R7</a><br>Description: 8" PRE-CRIMP 1858/19 RED      | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-V7</a><br>Description: 8" PRE-CRIMP 1858/19 VIOLET   | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-S9-D</a><br>Description: 8" PRE-CRIMP 1857/19 SLATE  | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-W9-D</a><br>Description: 8" PRE-CRIMP 1857/19 WHITE  | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-W7-D</a><br>Description: 8" PRE-CRIMP 1858/19 WHITE  | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-R7-D</a><br>Description: 8" PRE-CRIMP 1858/19 RED    | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-W7</a><br>Description: 8" PRE-CRIMP 1858/19 WHITE    | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-Y7</a><br>Description: 8" PRE-CRIMP 1858/19 YELLOW   | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-S9</a><br>Description: 8" PRE-CRIMP 1857/19 SLATE    | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-10-A7</a><br>Description: 10" PRE-CRIMP 1858/19 ORANGE  | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |
|  | Part#: <a href="#">1722533012-08-Y9</a><br>Description: 8" PRE-CRIMP 1857/19 YELLOW   | Manufacturers: <a href="#">Affinity Medical Technologies - a Molex company</a> | <div>RFQ</div> |

Related keywords for 1722533012-08-V9

|                                                                     |                                                      |                                                         |                                                    |
|---------------------------------------------------------------------|------------------------------------------------------|---------------------------------------------------------|----------------------------------------------------|
| Affinity Medical Technologies - a Molex company<br>1722533012-08-V9 | 1722533012-08-V9 distributor                         | 1722533012-08-V9 supplier                               | 1722533012-08-V9 price                             |
|                                                                     | 1722533012-08-V9 download datasheet                  | 1722533012-08-V9 datasheet                              | 1722533012-08-V9 stock                             |
| buy 1722533012-08-V9                                                | Affinity Medical Technologies - a Molex company      | Molex Connector Corporation 1722533012-08-V9            | Molex, LLC 1722533012-08-V9                        |
|                                                                     | 1722533012-08-V9                                     | Molex - Temp Flex 1722533012-08-V9                      | Oplink Communications, LLC 1722533012-08-V9        |
| Brad Harrison - a Molex company 1722533012-08-V9                    | Beau Interconnect - a Molex company 1722533012-08-V9 | Bittware - a Molex company 1722533012-08-V9             | FCT Electronics - a Molex company 1722533012-08-V9 |
|                                                                     | GWConnect - a Molex company 1722533012-08-V9         | Interconnect Systems - a Molex company 1722533012-08-V9 | Molex FCT Electronics 1722533012-08-V9             |
| Molex Brad Harrison 1722533012-08-V9                                | Molex Woodhead 1722533012-08-V9                      | Molex Flamar 1722533012-08-V9                           | Molex Phillips-Medisize 1722533012-08-V9           |
| Molex Interconnect Systems 1722533012-08-V9                         | Molex Affinity Medical Technologies 1722533012-08-V9 | Molex Beau Interconnect 1722533012-08-V9                | Molex PolyMicro Technologies 1722533012-08-V9      |
|                                                                     |                                                      | Molex NuCurrent 1722533012-08-V9                        | Molex SensorCon 1722533012-08-V9                   |
| Molex Triton Manufacturing 1722533012-08-V9                         | Molex Bittware 1722533012-08-V9                      | Molex Nallatech 1722533012-08-V9                        | Molex Nistica 1722533012-08-V9                     |

|                                                      |                                                           |                                                         |                                              |
|------------------------------------------------------|-----------------------------------------------------------|---------------------------------------------------------|----------------------------------------------|
| Molex 1722533012-08-V9                               | NuCurrent - a Molex company 1722533012-08-V9              | Nallatech - a Molex company 1722533012-08-V9            | Nistica - a Molex company 1722533012-08-V9   |
| Phillips-Medysize - a Molex company 1722533012-08-V9 | PolyMicro Technologies - a Molex company 1722533012-08-V9 | SensorCon - a Molex company 1722533012-08-V9            | Temp-Flex - a Molex company 1722533012-08-V9 |
|                                                      |                                                           | Triton Manufacturing - a Molex company 1722533012-08-V9 | Woodhead - a Molex company 1722533012-08-V9  |
|                                                      |                                                           |                                                         |                                              |

